



Colourless



F

Water- and oil-repellent



Highly wear-resistant




Achromatic





Reflection curve



Product structure



Quality tests

Properties

Optical		
	Reflection ≤3%	—
	UV filter	_
	Achromatic	

Mechanical

Hard	•
Highly wear-resistant	٠
Anti-static	٠

Chemical

Water- and oil-repellent	•
Resistant to solvents	٠
Low-temperature process	—
Alkali-resistant	—
Resistant to ultrasound	٠

Geometric

Drilling, opening	_	
Highly convex	٠	

Mechanical resistance	Description	Results
Vibratory finishing (chemical- mechanical polishing)	$2 \ hours in a vibrating bowl with cylindrical ceramic abrasives (Al_2O_3)$	Class 1
Scratch-Test	Test to establish the hardness and adhesiveness of the coating using a diamond tip and variable loads	manner internet tither the
Compliance with DIN-ISO 9211-4-02-02	Rapid pull off testing (~1 s) using adhesive tape 12-13 mm wide. Adhesive strength >9.8 N per 25 mm	Class 0 "OK"
Chemical resistance		
"Thermal shock" test DIN-ISO 9211-4-04-12	3×2-minute cycles in de-mineralised water at 100 °C. 1 minute in de-mineralised water at room temperature	Class 0 "OK"
Solubility DIN-ISO 9211-4-04-06	96 hours in a saline solution (45 g/l) at room temperature	Class 0 "OK"
Solvent resistance DIN-ISO 9211-3-12-3 and 9022-87-04-1	$1h$ hour in acetone (CH_3COCH_3) at room temperature	Class 0 "OK"

Din-150 9211-3-12-3 and alure OK 9022-87-04-1 Alkali bath (3% "Galvex") for 2 hours at 60 °C and Class 1 Resistance to ultrasound Alkali bath (3% "Galvex") for 2 hours at 60 °C and Class 1 "Minimal colour change" 7 days in an artificial sweat solution Class 0 DIN-ISO 9022-86-02-1 "OK "OK

esistance to chinatic conditions		
Cold	16 hours at -55 °C	Class O
DIN-ISO 9022-10-09-1		"OK"
Dry heat	6 hours at 85 °C with relative humidity of <40%	Class O
DIN-ISO 9022-11-06-1		"OK"

Cleaning instructions

Please observe the following instructions on how to clean glasses that have been treated with an anti-reflective coating:

 If cleaning by hand, we recommend wearing single-use, powder-free latex gloves and using a lint-free cloth.

 Glasses can be cleaned easily in a neutral solution using ultrasound and at temperatures of up to 60 °C. – Use one of the following approved cleaning alcohols: Ethanol with 5% isopropyl alcohol (F25-A+IPA), Alcosuisse, Bern Ethanol Absolute A 15 0 (02883), Sigma-Aldrich, Buchs SG

